

**Aprotinin Powder, Lyoph.**

**800277**

**1 Identification of the substance/mixture and of the company/undertaking**

**1.1 Product identifier**

Trade name Aprotinin Powder, Lyoph.  
Article number 800277

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Application of the substance / the mixture For research use only

**1.3 Details of the supplier of the safety data sheet**

Manufacturer/Supplier DSM Nutritional Products Ltd Branch Pentapharm  
Dornacherstrasse 112  
CH- 4147 Aesch BL / Switzerland  
Phone +41 61 706 48 48  
Fax +41 61 706 48 00  
Further information obtainable from pentapharm@pentapharm.com

**1.4 Emergency telephone number**

Tox Info Suisse : +41 44 251 51 51 (24h; 7 days a week)

**2 Hazards identification**

**2.1 Classification of the substance or mixture**

Classification according to Regulation (EC) No 1272/2008



Health hazard

Resp. Sens. 1 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.



Skins Sens. 1 H317 May cause an allergic skin reaction

**2.2 Label elements**

Labelling according to Regulation (EC) No 1272/2008

The substance is classified and labelled according to the CLP regulation.

Hazard pictograms



GHS08

# SAFETY DATA SHEET

According to 1907/2006/EC, Article 31



## Aprotinin Powder, Lyoph.

800277

### Signal word

Danger

### Hazard-determining components of labelling

Aprotinin

### Hazard statements

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H317 May cause an allergic skin reaction.

### Precautionary statements

P284 In case of inadequate ventilation wear respiratory protection.

P261 Avoid breathing dust.

P280 Wear protective gloves / protective clothing.

P342+P311 If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

P304+P341 IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

## 2.3 Other hazards

### Results of PBT and vPvB assessment

PBT: According to current knowledge, the product contains no PBT substance.



vPvB: According to current knowledge, the product contains no vPvB substance.

## 3 Composition/information on ingredients

### 3.2 Chemical characterization

#### Description

Mixture of substances listed below with nonhazardous additions.

Dangerous components		
CAS: 9087-70-1	Aprotinin	50 - 100%
EINECS: 232-994-9	 Resp. Sens. 1, H334;  Skin Sens. 1, H317	

#### SVHC

According to current knowledge no SVHC contained.

#### Additional information

For the wording of the listed risk phrases refer to section 16.

## 4 First aid measures

### 4.1 Description of first aid measures

#### General information

No special measures required

# SAFETY DATA SHEET

According to 1907/2006/EC, Article 31



## Aprotinin Powder, Lyoph.

800277

### After inhalation

If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

### After skin contact

Wash with plenty of soap and water.

If skin irritation or rash occurs: Get medical advice/attention.

### After eye contact

Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

If eye irritation persists: Get medical advice/attention.

### After swallowing

Immediately call a POISON CENTRE/doctor.

## 4.2 Most important symptoms and effects, both acute and delayed

May cause an allergic skin reaction.

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

## 4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

## 5 Firefighting measures

### 5.1 Extinguishing media

Suitable extinguishing media

CO<sub>2</sub>, foam, dry powder, or water.

### 5.2 Specific hazards arising from the substance or mixture

No further relevant information available.

### 5.3 Advice for firefighters

Protective equipment

Wear self-contained respiratory protective device.

Wear protective clothing to prevent contact with skin and eyes.

## 6 Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Avoid skin and eye contact and particle inhalation.

### 6.2 Environmental precautions

Do not allow product to reach sewage system or any water course.

### 6.3 Methods and material for containment and cleaning up

Collect in suitable containers.

Version 01 - Revision 27.09.2016  
800277\_SDS\_EN

# SAFETY DATA SHEET

According to 1907/2006/EC, Article 31



## Aprotinin Powder, Lyoph.

800277

Clean the affected area carefully and ensure adequate ventilation.

Avoid dust information.

### 6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

## 7 Handling and storage

### 7.1 Precautions for safe handling

Avoid skin and eye contact particle inhalation

Information about fire - and explosion protection

No special measures required.

### 7.2 Conditions for safe storage, including any incompatibilities

**Storage**

**Requirements to be met by storerooms and receptacles**

No special requirements.

**Information about storage in one common storage facility**

Not required.

**Further information about storage conditions**

Store the product protected from moisture and light between +2 to +8 °C.

### 7.3 Specific end use(s)

No further relevant information available.

## 8 Exposure controls/personal protection

### 8.1 Control parameters

**Ingredients with limit values that require monitoring at the workplace**

A limit value has not been defined.

**Additional information**

The lists valid during the making were used as basis.

The professional use of this product by young people is restricted or prohibited altogether. The corresponding legal bases and precise regulations are set out in section 15.

### 8.2 Exposure controls

**Personal protective equipment**

**General protective and hygienic measures**

Ensure good ventilation/exhaustion at the workplace.

Avoid contact with the eyes and skin.

Keep away from foodstuffs, beverages and feed.

Immediately remove all soiled and contaminated clothing

**Aprotinin Powder, Lyoph.**

**800277**

Wash hands before breaks and at the end of work.

**Respiratory protection**

In case of inadequate ventilation wear respiratory protection.

Respiratory protection must be compliant with the specifications of Directive 89/686/EEC and be in accordance with the EN standards resulting thereof.

**Protection of hands**



Protective gloves

**Material of gloves**

Protective gloves must be compliant with the specifications of EU Directive 89/686/EEC and be in accordance with the standard EN374. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation.

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

**Eye protection**



Safety glasses

Protective goggles must be compliant with the specifications of Directive 89/686 / EEC and be in accordance with the standard EN166.

**Limitation and supervision of exposure into the environment**

Do not allow product to reach sewage system or any water course.

**9 Physical and chemical properties**

**9.1 Information on basic physical and chemical properties**

General Information	
<b>Appearance</b>	
Form	Voluminous lyophilisate
Colour	Almost white
Odour	Odourless
Odour threshold	Not determined.
pH-value	Not determined.
<b>Change in condition</b>	
Melting point/Melting range	Not determined.

# SAFETY DATA SHEET

According to 1907/2006/EC, Article 31



## Aprotinin Powder, Lyoph.

800277

Boiling point/Boiling range	Not determined.
Flash point	Not applicable.
Flammability (solid, gaseous)	Not determined.
Ignition temperature	Not determined.
Decomposition temperature	Not determined.
Self-igniting	Product is not self-igniting.
Danger of explosion	Not determined
Explosion limits	
Lower	Not determined.
Upper	Not determined.
Oxidising properties	Not determined.
Vapour pressure	Not applicable.
Density	
Relative density at 20 °C	Not determined.
Vapour density	Not determined.
Evaporation rate	Not applicable.
Solubility in / Miscibility with	
Water	Soluble.
Partition coefficient (n-octanol/water)	Not determined.
Viscosity	
Dynamic	Not applicable.
Kinematic	Not applicable.
Solvent content	
VOC (EC)	Not determined.
VOCV (CH)	Not determined.

### 9.2 Other information

No further relevant information available.

## 10 Stability and reactivity

### 10.1 Reactivity

No further relevant information available.

### 10.2 Chemical stability

No information available.

**Thermal decomposition / conditions to be avoided**

No decomposition if used according to specifications.

### 10.3 Possibility of hazardous reactions

No dangerous reactions known.

## **Aprotinin Powder, Lyoph.**

**800277**

### **10.4 Conditions to avoid**

Heat may cause this product to decompose.

### **10.5 Incompatible materials**

None known

### **10.6 Hazardous decomposition products**

Nature of decomposition products not known.

## **11 Toxicological information**

### **11.1 Information on toxicological effects**

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

#### **Skin corrosion/irritation**

Based on available data, the classification criteria are not met.

#### **Serious eye damage/irritation**

Based on available data, the classification criteria are not met.

#### **Respiratory or skin sensitisation**

May cause an allergic skin reaction.

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

#### **Germ cell mutagenicity**

Based on available data, the classification criteria are not met.

#### **Carcinogenicity**

Based on available data, the classification criteria are not met.

#### **Reproductive toxicity**

Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

## **12 Ecological information**

### **12.1 Toxicity**

#### **Aquatic toxicity**

No further relevant information available.

### **12.2 Persistence and degradability**

No further relevant information available.

According to 1907/2006/EC, Article 31



## 800277

No further relevant information available.

No further relevant information available.

PBT

According to current knowledge, the product contains no PBT substance.

vPvB

According to current knowledge, the product contains no vPvB substance.

No further relevant information available.

### 13.1 Waste treatment methods

Recommendation: Disposal must be made according to official regulations.

European waste catalogue

The Swiss waste code is identical.

In Switzerland codes for hazardous waste are marked with an “S” (“Sonderabfall”) instead of an asterisk.

18 00 00	WASTES FROM HUMAN OR ANIMAL HEALTH CARE AND/OR RELATED RESEARCH (except kitchen and restaurant wastes not arising from immediate health care)
18 01 00	wastes from natal care, diagnosis, treatment or prevention of disease in humans
18 01 06*	chemicals consisting of or containing dangerous substances

\* = hazardous waste

### Uncleaned packaging

Recommendation: Disposal must be made according to official regulations.

Recommended cleansing agents: Water, if necessary together with cleansing agents.

### 14.1 UN-number

ADR, ADN, ADN, IMDG, IATA Void

ADR, ADN, ADN, IMDG, IATA

Void

ADR, ADN, ADN, IMDG, IATA

Class	Void
-------	------



According to 1907/2006/EC, Article 31



## 800277

ADR, ADN, IMDG, IATA

Void

## Marine pollutant

No

Not applicable.

Not applicable.

### Transport/Additional information

ADR/ADN

Remarks

Not classified.

IMDG

Remarks

Not classified.

IATA

Remarks

Not classified.

## UN "Model Regulation"

Void

## 15 Regulatory information

## Directive 2012/18/EU

## Named dangerous substances - ANNEX I

None of the ingredients is listed.

## National regulations

**Information about limitation of use:**

Article 4 para. 4 of the Ordinance on the protection of young people in the workplace (SR 822.115) and Article 1lit. f of the EAER regulation on hazardous work and young people (SR 822.115.2) Young people undergoing basic vocational training may only work with this product if the relevant training ordinance makes provision for them to do so with a view to enabling them to achieve their training objectives and if the preconditions for the training plan have been met and the applicable age restrictions have been complied with. Young people who are not completing any basic vocational training are not permitted to work with this product. Employees of either sex who are under 18 years old are classed as young people.

### Waterhazard class

Water hazard class 1 (Self-assessment): slightly hazardous to water.

A Chemical Safety Assessment has not been carried out.

# SAFETY DATA SHEET

According to 1907/2006/EC, Article 31



**Aprotinin Powder, Lyoph.**

**800277**

## 16 Other information

All information and instructions provided in this Safety Data Sheet are based on the current state of scientific and technical knowledge at the date indicated on this Safety Data Sheet. DSM Nutritional Products Ltd Branch Pentapharm shall not be held responsible for any defect in the product covered by this Safety Data Sheet, should the existence of such a defect not be detectable considering the current state of scientific and technical knowledge.

### Relevant phrases

H317 May cause an allergic skin reaction.

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

### Abbreviations and acronyms

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International

Transport of Dangerous Goods by Rail)

ICAO: International Civil Aviation Organisation

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of

Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOCV: Lenkungsabgabe auf flüchtigen organischen Verbindungen, Schweiz (Swiss Ordinance on volatile organic compounds)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Resp. Sens. 1: Sensitisation - Respirat., Hazard Category 1

Skin Sens. 1: Sensitisation - Skin, Hazard Category 1